



NaviAid™ AB

YOUR GI ENDOSCOPE CAN
ACCOMPLISH SO MUCH MORE



- ⦿ Deep ileoscopy / upper enteroscopy
- ⦿ Fast and safe advancement in the small bowel
- ⦿ Introduced upon need, through the endoscope's instrument channel
- ⦿ Enables Controlled Withdrawal™ of the endoscope

ENDOSCOPY, THE SMART WAY

smart[⦿]
medical systems



NaviAid™ AB

YOUR GI ENDOSCOPE CAN ACCOMPLISH SO MUCH MORE

The NaviAid™ AB (Advancing Balloon) device enables you to carry out fast and easy ileoscopy and upper enteroscopy procedures with your standard endoscope. Used on demand, no pre-mounting or preparation prior to the procedure is necessary. The device is introduced through the tool channel of the endoscope, on the physician's decision, and advanced ahead of the endoscope. Once inflated in its position

preceding the endoscope and anchored in the bowel, the device becomes a "rail" on which the endoscope can be advanced, enabling rapid, safe intubation in the small bowel. Multiple advancement steps may be performed as needed, and the device may be removed and replaced at will, permitting biopsies and therapy to be conducted.

ADDITIONAL KEY ADVANTAGES:

- Works with your standard endoscope
- Maintains the endoscope's biopsy and therapy capabilities
- Intuitive, easy to operate

NaviAid™ AB – PRODUCT SPECIFICATIONS:

Catheter Length	3,500mm
Balloon Diameter	40mm
Inflated Balloon Set Pressure	60mbar
Set Pressure Tolerance	±10mbar
Balloon Inflation System- Power Supply	100-240VAC; 50-60Hz



ABOUT SMART MEDICAL

Your professional partner in a variety of endoscopy applications, SMART is a pioneer in the development and manufacture of world-leading medical devices. SMART's NaviAid™ product family enhances the capabilities and range of standard endoscopes in an intuitive, cost-effective manner. SMART is headquartered in Israel's hi-tech region, and operates in the United States via its wholly-owned subsidiary, SMART GI Inc. SMART is ISO 13485 certified, and its globally-distributed commercial products have gained FDA clearance (in the USA) and the CE Mark (in Europe).



ENDOSCOPY, THE SMART WAY

Headquarters: 5 HaNofar st, P.O Box 4123 Raanana, 4366404 Israel. Tel:+972-9-7444321, Fax:+972-9-7444543, info@smartmedsys.com • US Offices: SMART GI Inc., 540 Frontage Road, Suite 2265, Northfield, IL 60093, Tel: 1-847-920-5442, Fax: 1-800-230-5823 us@smartmedsys.com

www.smartmedsys.com

© 2019 Smart Medical Systems Ltd. All rights reserved.