



G-EYE® Colonoscope

YOUR STANDARD COLONOSCOPE
CAN ACCOMPLISH SO MUCH MORE



- Compatible with selected 510(k) cleared brands and models
- Reusable reprocessable balloon
- Maintaining colonoscope's field of view
- Endoscope anchoring during therapy

ENDOSCOPY, THE SMART WAY

smart^o
medical systems

The G-EYE® Colonoscope is constructed of a standard, 510(k) cleared colonoscope, which is remanufactured to include a reusable and reprocessable balloon at its bending section. During colonoscopy, withdrawal of the colonoscope with the G-EYE® balloon moderately inflated flattens colon folds, centralizes the optics within the colonic lumen, and assists in controlling the view-field and positioning of the colonoscope.

Upon need, the G-EYE® balloon can be inflated at anchoring pressure to stabilize the colonoscope tip during endoscopic intervention and therapy.

The G-EYE® Colonoscope further assists in avoiding colonoscope slippage in the bowel, thus ensuring a steady and Controlled Withdrawal™. The G-EYE® Colonoscope does not interfere with the use of endoscopic accessories that may be introduced through the colonoscope's instrument channel when necessary. The inflation of the G-EYE® balloon is controlled by the NaviAid™ SPARKC inflation system.

G-EYE® - PRODUCT SPECIFICATIONS:

Available endoscopes	Selected 510(k) cleared brands and models
Balloon Diameter	Up to 60mm
Inflated Balloon Set Pressure	70mbar
Controlled Withdrawal™ pressure	3 intermediate pressure levels
SPARKC Power Supply	100V - 240VAC; 50-60Hz



ABOUT SMART MEDICAL

SMART Medical Systems Ltd. is a pioneer in the development and manufacturing of medical devices, specializing in the field of gastroenterology and gastrointestinal (GI) endoscopy. SMART's unique approach is to use available brand name endoscopes and address key challenges in contemporary endoscopy. SMART Medical is headquartered in Israel's hi-tech region, and operates in the United States via its wholly-owned subsidiary, SMART GI Inc. SMART Medical is an ISO 13485 certified company.



ENDOSCOPY, THE SMART WAY

Headquarters: 5 HaNofar st, P.O Box 4123 Raanana, 4366404 Israel. Tel:+972-9-7444321, Fax:+972-9-7444543, info@smartmedsys.com • US Offices: SMART GI Inc., 540 Frontage Road, Suite 2265, Northfield, IL 60093, Tel: 1-847-920-5442, Fax: 1-800-230-5823 us@smartmedsys.com